

The Mountain Pine Beetle (*Dendroctonus ponderosae*)

Synergy Semiochemicals Lure (#3125) with Terpinolene

The Synergy Lure (#3125) is made up of 3 components:

- A) Exo brevicomin Flexlure
→comes in a brown bag with C
- B) Terpinolene Polybottle
→comes in a white bag
- C) Trans verbenol bubble cap



Field Life: Contents slowly release through the plastic pouches and last approximately **90 days @ 20° C.**

Storage: Stored un-used lures in a freezer.

Lure Hanging: Thread twist-tie (provided) through the hole-punch near the top of the bags. Tether components together to the outside of the trap. We recommend tethering the lures toward the bottom third of the trap.

Important Note! Do NOT open pouches or puncture any of the plastic bubbles. These pouches are part of the controlled release system and protect the lure's active ingredients from the outside elements.

Safety and Disposal: Components are generally non-toxic chemicals with little risk in handling them. However, users should use gloves or wash their hands after handling and keep lures out of reach from children and pets. Spent lures may be disposed of in normal trash.

General Info: This is a funnel lure used to attract the beetle. Traps are used as part of a monitoring or Integrated Pest Management strategy.

Repellents and other products: If you want to **repel** Mountain Pine Beetles please refer to our Synergy Shield Verbenone Pouches (#3413 and #3414).

Synergy Semiochemicals Corp.

1-604- 454-1122 | www.semiochemical.com | synergy@semiochemical.com

