

The Asian Longhorn Beetle (*Anoplophora glabripennis*) Synergy Semiochemicals Lure (#3432)

The Asian Longhorn Beetle (ALB)
Synergy Semiochemicals Lure (#3432)
is made up of 3 components:

- A) ALB aldehyde
→comes in clear mini-bubble
- B) ALB alcohol
→comes in a white bubble
- C) ALB kairomone blend
→comes in blue bag



Important Note! Both the ALB aldehyde and alcohol must be switched out and replaced by new components after 4-6 weeks (30-40 days) after deployment. Extra ALB aldehyde (#3211) and ALB alcohol (#3240) bubbles must be ordered in addition to the ALB Synergy Combo Lure

Field Life: Pouch contents slowly release through the plastic and last approximately **90 days** @ 20° C. Bubble caps release within 45 days.

Storage: Stored un-used lures in a freezer.

Lure Hanging: Thread twist-tie (provided) through the hole-punch near the top of the bags. Tether components together to the outside of the trap. We recommend tethering the lures toward the bottom third of the trap.

Important Note! Do NOT open pouches or puncture any of the plastic bubbles. These pouches are part of the controlled release system and protect the lure's active ingredients from the outside elements.

Safety and Disposal: Components are generally non-toxic chemicals with little risk in handling them. However, users should use gloves or wash their hands after handling and keep lures out of reach from children and pets. Spent lures may be disposed of in normal trash.

Synergy Semiochemicals Corp.

1-604- 454-1122 | www.semiochemical.com | synergy@semiochemical.com

